

PathoSEEK[®] Total Aerobic Bacteria Count Assay v2 with SenSATIVAx[®] DNA Purification Protocol for Detection in Cannabis Flower and MIP Matrices

Method Developer Validation

Table of Contents:

Abstract	3
Materials	4
Inclusivity and Exclusivity	5
Table 1: Inclusivity Results, PathoSEEK® Total Aerobic Count Assay	5
Table 2: Exclusivity Results, PathoSEEK® Total Aerobic Count Assay	7
Cq to CFU Equation for Flower Samples	8
Proficiency Testing / Certified Reference Materials	8
Table 3: CRM Results Flower, PathoSEEK® Total Aerobic Count Assay	9
Table 4: Detection of 100 cfu in Oil, PathoSEEK® Total Aerobic Count Assay	9
Table 5: CRM Results MIP, PathoSEEK® Total Aerobic Count Assay	10
Conclusions	10

Abstract

Background:

Total Aerobic Bacteria encompasses many species. Bacteria can cause deterioration and decomposition of cannabis, and certain species of aerobic bacteria, such as Shiga Toxin producing *E. coli*, can cause infections in humans. Current regulations allow cannabis flower and cannabis products to contain a limit of total aerobic bacteria. The PathoSEEK® Total Aerobic Bacteria Count, or Total Aerobic Count (TAC), Assay with SenSATIVAx® DNA Purification Protocol is designed to detect all aerobic bacteria in a single qPCR (Quantitative Polymerase Chain Reaction) in cannabis flower, hemp flower, cannabis concentrates, infused edibles, and infused non-edibles.

Objective:

To evaluate the PathoSEEK® Total Aerobic Count Detection Assay v2, using the SenSATIVAx® Flower and MIP DNA purification protocols for the enumeration of total aerobic bacteria in cannabis flower (delta 9-tetrahydrocannabinol >0.3%; 1g), and for presence/absence detection in marijuana-infused products (MIP).

Results:

The inclusivity and exclusivity results demonstrated the high specificity of the PathoSEEK® Total Aerobic Count method in differentiating target organisms, prevalent in cannabis flower and infused products, from non-target organisms. The SenSATIVAx® flower DNA Purification kit and PathoSEEK® Total Aerobic Count Detection Assay v2 was validated by constructing an enumeration curve and conversion equation using eleven distinct bacterial species and plating on 3M™ Petrifilm™ RAC (Rapid Aerobic Count) Plates. Subsequent to the generation of this curve, a Certified Reference Materials (CRMs) from NSI, Quantitative APC in Hemp and EB in Hemp, were procured and analyzed utilizing the SenSATIVAx flower DNA Purification kit and PathoSEEK Total Aerobic Count qPCR Detection Assay v2. Resulting Cq values were converted to CFU/g utilizing the conversion equation and were compared to results obtained via Neogen Petrifilm™ RAC Petrifilms and to the NSI value provided on the CRM Certificate of Analysis. Results generated by qPCR were comparable to the RAC Petrifilm results and aligned with the specifications presented by NSI.

Materials

SenSATIVAx® Flower & Leaf DNA Purification Kit Components - P/N 420001

Component Name	Qty Provided	Storage Conditions
MGC Cell Lysis Buffer	1 Bottle (12 mL)	RT (20–28°C)
MGC Binding Buffer	1 Bottle (48 mL)	Refrigerate (2-8 °C)
MGC Elution Buffer	1 Bottle (12 mL)	RT (20–28°C)

SenSATIVAx® Infused Product DNA Purification Kit Components - P/N 420004

Component Name	Qty Provided	Storage Conditions
SenSATIVAx® Solution A	1 Bottle (350 mL)	RT (20–28°C)
SenSATIVAx® Solution B	1 Bottle (25 mL)	RT (20–28°C)
MGC Binding Buffer	1 Bottle (48 mL)	Refrigerate (2-8 °C)
MGC Elution Buffer	1 Bottle (12 mL)	RT (20–28°C)

PathoSEEK® Internal Control - P/N 420337

Component Name	Qty Provided	Storage Conditions
Internal Control	1 Tube (50 µL)	-15 to -20 °C

PathoSEEK® Total Aerobic Bacteria Count (TAC) Detection Assay v3 Kit - P/N 420541

Component Name	Qty Provided	Storage Conditions
PathoSEEK® Total Aerobic Count Detection Assay v2	1 Tube (200 µL)	-15 to -20 °C
qPCR Master Mix Kit v3 - Reaction Buffer	1 Tube (160 µL)	-15 to -20 °C
qPCR Master Mix Kit v3 - Nuclease Free Water	2 Tubes (1 mL)	-15 to -20 °C
qPCR Master Mix Kit v3 - Master Mix v3	1 Tube (750 µL)	-15 to -20 °C

Optional: Grim Reef Free DNA Removal Kit - P/N 420145

Component Name	Qty Provided	Storage Conditions
GR Enzyme	1 Bottle (2.5 mL)	-15 to -20 °C
GR Buffer	1 Bottle (12.5 mL)	-15 to -20 °C

Note: Actual fill volumes include overage

Inclusivity and Exclusivity

Methodology

For the inclusivity evaluation, 30 bacterial strains were assessed. Target strains were either cultured in Tryptic Soy Broth for 24 hours at 37°C, followed by DNA extraction, or purified DNA from ATCC was utilized. For exclusivity, 30 organisms were evaluated. Target strains were either cultured under optimal growth conditions for the organism, followed by DNA extraction, or purified DNA from ATCC was utilized. Inclusivity and exclusivity cultures were randomized, blind-coded, and analyzed using the PathoSEEK® Total Aerobic Count method.

Results

Of the 30 inclusivity strains tested, 30 were correctly detected by the method. Of the 30 exclusivity strains tested, all 30 were correctly excluded. Tables 1 and 2 present a summary of the results.

Table 1: Inclusivity Results, PathoSEEK® Total Aerobic Count Assay

Species	ATCC#	PathoSEEK®
		TAC Result
<i>Acinetobacter baumannii</i>	19606	Detected
<i>Aeromonas hydrophila</i>	7966	Detected
<i>Burkholderia multivorans</i>	17616	Detected
<i>Bacillus subtilis</i>	11774	Detected
<i>Citrobacter braakii</i>	3037	Detected
<i>Citrobacter koseri</i>	25408	Detected
<i>Edwardsiella tarda</i>	23672	Detected

<i>Enterobacter aerogenes</i>	13048	Detected
<i>Enterobacter cloacae</i>	13047	Detected
<i>Erwinia rhapontici</i>	29290	Detected
<i>Escherichia coli</i>	25922	Detected
<i>Escherichia coli O157:H7</i>	35150	Detected
<i>Escherichia hermannii</i>	700368	Detected
<i>Escherichia vulneris</i>	33821	Detected
<i>Hafnia alvei</i>	51873	Detected
<i>Klebsiella oxytoca</i>	51983	Detected
<i>Klebsiella pneumonia</i>	BAA-2146	Detected
<i>Listeria monocytogenes</i>	7647	Detected
<i>Morganella morganii</i>	25829	Detected
<i>Pantoea agglomerans</i>	43348	Detected
<i>Proteus mirabilis</i>	43071	Detected
<i>Pseudomonas aeruginosa</i>	15442 & 35552	Detected
<i>Pseudomonas fluorescens</i>	13525	Detected
<i>Pseudomonas putida</i>	47054	Detected
<i>Ralstonia insidiosa</i>	49129	Detected
<i>Rahnella species</i>	33991	Detected
<i>Salmonella enterica</i>	13311	Detected
<i>Stenotrophomonas maltophilia</i>	13637	Detected
<i>Staphylococcus aureus</i>	12600	Detected
<i>Serratia marcescens</i>	27137	Detected

Table 2: Exclusivity Results, PathoSEEK® Total Aerobic Count Assay

Species	ATCC#	PathoSEEK®
		TAC Result
<i>Alternaria alternata</i>	6663	Not Detected
<i>Aspergillus alabamensis</i>	MYA-3633	Not Detected
<i>Aspergillus brasiliensis</i>	9642	Not Detected
<i>Aspergillus carbonarius</i>	MYA-4641	Not Detected
<i>Aspergillus caesiellus</i>	42693	Not Detected
<i>Aspergillus carneus</i>	13549	Not Detected
<i>Aspergillus clavatus</i>	1007	Not Detected
<i>Aspergillus deflectus</i>	62502	Not Detected
<i>Aspergillus flavus</i>	9643	Not Detected
<i>Aspergillus fumigatus</i>	1022	Not Detected
<i>Aspergillus japonicus</i>	16873	Not Detected
<i>Aspergillus nidulans</i>	38163	Not Detected

<i>Aspergillus niger</i>	13496	Not Detected
<i>Aspergillus oryzae</i>	42149	Not Detected
<i>Aspergillus parasiticus</i>	56775	Not Detected
<i>Aspergillus pseudoterreus</i>	10020	Not Detected
<i>Aspergillus terreus</i>	1012	Not Detected
<i>Aspergillus tubigenensis</i>	1004	Not Detected
<i>Aspergillus ustus</i>	1041	Not Detected
<i>Aspergillus versicolor</i>	11730	Not Detected
<i>Candida albicans</i>	10231	Not Detected
<i>Cryptococcus laurentii</i>	18803	Not Detected
<i>Cryptococcus neoformans</i>	208821	Not Detected
<i>Fusarium oxysporum</i>	62506	Not Detected
<i>Fusarium solani</i>	52628	Not Detected
<i>Mucor luteus</i>	28932	Not Detected
<i>Penicillium chrysogenum</i>	18476	Not Detected
<i>Penicillium rubens</i>	11709	Not Detected
<i>Rhizopus stolonifera</i>	14037	Not Detected
<i>Talaromyces pinophilus</i>	11797	Not Detected

Cq to CFU Equation for Flower Samples

1. The equation utilized for quantifying colony-forming units per gram (CFU/g) was established through quantitative polymerase chain reaction (qPCR) analysis of eleven bacterial organisms, followed by their subsequent cultivation on Petrifilm Coliform Count (CC) plates. Both qPCR and plating protocols were executed in triplicate across a series of six serial dilutions. Mean values were computed and depicted on a scatter graph, with qPCR cycle quantification (Cq) data presented on the x-axis and the base-10 logarithm of plating data on the y-axis. A linear regression analysis was performed to determine the best-fit line, resulting in the equation: $y = -0.2383x + 10.005$, where y is the log₁₀ CFU/g and x is the Cq value.
2. Utilize the following linear equation to convert Cq (x) values to Log CFU (y): $y = -0.2383x + 10.005$.
3. Perform an inverse logarithmic transformation of Y to obtain CFU/g.
4. Multiply the derived CFU/g value by the sample's upfront dilution factor (x 20) in TSB to determine the final CFU/g.

Empirical validation has confirmed that this derived equation produces comparable results across the

Agilent AriaMX, BioRad CFX96, and BioMolecular Systems MIC quantitative PCR instruments.

Proficiency Testing / Certified Reference Materials

Methodology - Flower Matrix

Two CRMs from NSI, APC in Hemp (catalog # FM-725) and EB in Hemp (catalog # FM-730), were analyzed following the PathoSEEK Total Aerobic Count User Guide v3. qPCR was performed in triplicate on an Agilent Aria, Bio-Rad CFX96, and BioMolecular Systems MIC. Cq values were converted into CFU/g values and compared with the expected range of the CRMs provided by NSI.

Table 3: CRM Results Flower, PathoSEEK® Total Aerobic Count Assay

Content	Instrument	Average Cq Fam	CFU
APC in Hemp	Agilent Aria	27.92	45,029
APC in Hemp	Biorad CFX96	28.62	30,565
APC in Hemp	BMS MIC	28.19	38,730
EB in Hemp	Agilent Aria	27.68	51,273
EB in Hemp	Biorad CFX96	28.24	37,764
EB in Hemp	BMS MIC	27.88	45,826

Methodology, Non-Flower Matrices

To ascertain our ability to detect Total Aerobic Bacteria in MIPs, we acquired a *S. arizonae* EZ Accusshot from Microbiologics (Ref # 0901A). Following the provided instructions, we inoculated 1g of Oil with 100 cfu of *S. arizonae*, 2.4 mL TSB and enriched the sample at 37°C for 16 hours. Post enrichment, three samples were processed and analyzed on all three qPCR instruments to determine if this concentration of Total Aerobic Bacteria could be detected.

Results

Each instrument detected the presence of Total Aerobic Bacterial DNA.

Table 4: Detection of 100 cfu in Oil, PathoSEEK® Total Aerobic Count Assay

Content	Instrument	Detected
S. arizonae in Oil	Agilent Aria	Yes
S. arizonae in Oil	Biorad CFX96	Yes
S. arizonae in Oil	BMS MIC	Yes

Methodology, Non-Flower Matrices Continued

To further test the assays ability to detect aerobic bacteria in non flower matrices, two CRMs from NSI, S. aureus in Hemp Oil (catalog # FM-832) and STEC in Edible (catalog # CMQC-029), were analyzed following the PathoSEEK Total Aerobic Count User Guide v3. Samples were enriched in TSB for 16 hours at 37°C followed by DNA purification and qPCR. qPCR was performed in triplicate on an Agilent Aria, Bio-Rad CFX96, and BioMolecular Systems MIC. A Cq value indicates a positive result and the sample material should be plated for enumeration.

Results

Each instrument called within the expected range of each CRM.

Table 5: CRM Results MIP, PathoSEEK® Total Aerobic Count Assay

Content	Instrument	Detected
S. aureus in Hemp oil	Agilent Aria	Yes
S. aureus in Hemp oil	Biorad CFX96	Yes
S. aureus in Hemp oil	BMS MIC	Yes
STEC in Edible	Agilent Aria	Yes
STEC in Edible	Biorad CFX96	Yes
STEC in Edible	BMS MIC	Yes

Results

Each instrument detected the presence of Total Aerobic Bacterial DNA.

Conclusions

The PathoSEEK® Total Aerobic Count qPCR Assay v2 with SenSATIVAx® DNA Purification is a rapid, alternative method to traditional plating procedures for the detection of aerobic bacteria in cannabis flower and cannabis infused products. The method produced comparable results to 3M Petrifilm EB and CC plates for the enumeration of aerobic bacteria in cannabis flower.

REVISION HISTORY

Version	Date	Description
v1	September, 2021	<ol style="list-style-type: none">1. Update to User Guide Format2. Updated sample to media ratio used for homogenization of flower3. Update to qPCR Master Mix v34. Update to conversion equations
v2	November 2022	<ol style="list-style-type: none">1. Updated sample to media ratio used for homogenization of flower2. Update to conversion equation for flower matrix3. Update MIP detection to a presence absence, removal of conversion equation for MIP matrices.
v3	August 2025	<ol style="list-style-type: none">1. New kitted packaging format2. New assay version (v2) to include use with BMS Mic and Myra3. BMS Mic and Myra data analysis section4. Update to Internal Control from SCCG to IC

DISCLAIMER

This test was developed, and its performance characteristics determined by Medicinal Genomics Company, for laboratory use. Any deviations from this protocol are not supported by MGC.

This test has not been validated on remediated (irradiated, ozone treated, acid treated, hydrogen peroxide treated, etc.) samples. Samples that have undergone remediation may cause discordant results between plating methods and PathoSEEK® methods. When remediated samples produce a result above the action limit on qPCR, we recommend confirming viability with an approved plating method.

Results may vary based on laboratory conditions. Altitude and humidity are factors known to affect the growth of bacterial and fungal species.

LIMITED USE LABEL LICENSE

This product is covered by at least one or more claims of US patent applications, which are exclusively licensed to Medicinal Genomics Corporation. This product is sold strictly for the use of the buyer, and the buyer is not authorized to transfer this product [or any materials made using this product] to any third party.

© 2025 Medicinal Genomics Corporation. All rights reserved.

* All Trademarks are property of their respective owners.